



In-silico bifurcation stenting: plaque composition / bifurcation angle vs plaque shift

J.J. Wentzel
C. Chiastra, F. Iannaccone


Erasmus MC, Rotterdam, The Netherlands
Politecnico di Milano, Milano, Italy
Gent University, Gent, Belgium



European Research Council
Established by the European Commission

To study the influence of distal angle / plaque composition on plaque shift / carina shift because of provisional stenting

- ⇒ Create model with plaque included
- ⇒ 3 bifurcation geometries with different distal angles* are investigated
- ⇒ 2 different type of plaques:
 - fibrotic
 - lipid




Novel bifurcation phantoms for validation of quantitative coronary angiography algorithms. Girasis C, ... Serruys, Wentzel JJ. CCI, 2011

* **Distal angle** = $57.3^\circ \pm 10.0^\circ$ calculated on LAD, RCA, LCX (mainly LAD, 92.2%) (n = 153 patients) by Elsaban et al. 2013
Elsaban et al. J Invasive Cardiol 2013; 25:118-122

Bifurcation phantom studies

Designed bifurcation phantom




1. Diameter ratio : finet law
2. Length lesion
3. Location MLD
4. Bifurcation angle

Novel bifurcation phantoms for validation of quantitative coronary angiography algorithms. Girasis C, ... Serruys ,Wentzel JJ. CCI, 2011

Model : LAD / D1 Bifurcation

- Intima
- Media
- Adventitia
- Plaque



Arterial wall thickness (Tzucu 2001, Holzapfel 2005):

intima + media thickness = 0.1 lumen Diameter (healthy)


intima thickness/wall thickness = 0.26

media thickness/total wall thickness = 0.35


adventitia thickness/wall thickness = 0.39

Plaque burden > 40% at location of lumen narrowing (Glagov 1987)

Model : LAD / D1 Bifurcation



PB = Plaque burden (up to the external elastic membrane)




Computer simulations

Mechanical properties of plaque components


Arterial wall

Isotropic hyperelastic behavior with ideal plasticity to mimic vessel damage



Fibrotic / Lipid plaque

Isotropic hyperelastic behavior with ideal plasticity to mimic plaque rupture



Methods :Computer simulations



Stent and balloon

- **Multi-Link 8** (Abbott Laboratories, Abbott Park, IL, USA)

- Bare-metal stent, Co-Cr alloy
- Size: 3x18 mm


- **Balloon:**

- Modeled as a straight tube using Kiousis approach*
- Calibrated using the manufacturer compliance chart from 10 atm to 14 atm (nominal pressure = 10 atm, burst pressure = 18 atm)




Computer simulated provisional stenting

1. Insertion of 3x18mm Multi-Link 8 stent




2. Balloon expansion at 14 atm

(Final diameter = 3.22 mm)




3. Final configuration after balloon deflation



Note: Post-dilatation with a 3.5x12 balloon is needed (work in progress)

Results – fibrotic plaque

- Geometrical configuration pre/post






45 deg




57.5 deg

70 deg

PRE



POST





- Distal angle decreases after stent deployment
- Lumen shape of side branch is effected

Results: plaque shift

- Volumetric analysis

FIBROTIC PLAQUE




SB ostial compromise = lumen volume decrease in the SB ostial segment after MB stenting
Plaque shift = plaque volume increase in the SB ostial segment after MB stenting
Carina shift = vessel volume decrease in the SB ostial segment after MB stenting

- SB lumen : Distal angle ↑↑ ⇒ SB ostial compromise ↓↓
- SB plaque : Distal angle ↑↑ ⇒ Plaque shift ↓↓ (absent with distal angle = 70°)

Distal angle ↑↑ ⇒ Area Polygon of confluence ↑↑ ⇒ Plaque shift ↓↓

Results – lipid plaque

- Geometrical configuration pre/post






45 deg




57.5 deg

70 deg

PRE





POST



- Distal angle decreases after stent deployment
- Lumen shape of side branch is effected
- The decrease of distal angle is greater in the lipid cases than the fibrotic cases

Results: plaque shift

- Volumetric analysis





SB ostial compromise = lumen volume decrease in the SB ostial segment after MB stenting
Plaque shift = plaque volume increase in the SB ostial segment after MB stenting
Carina shift = vessel volume decrease in the SB ostial segment after MB stenting

- SB lumen : SB ostial compromise and plaque shift more severe in **lipid** cases
- SB plaque : Carina shift present only in the **lipid** case with **distal angle = 70°**
- **Lumen size of side branches marginally changes after stenting but lumen shape is effected**


Results: plaque shift

- Volumetric analysis





SB ostial compromise = lumen volume decrease in the SB ostial segment after MB stenting
Plaque shift = plaque volume increase in the SB ostial segment after MB stenting
Carina shift = vessel volume decrease in the SB ostial segment after MB stenting

- SB lumen : SB ostial compromise and plaque shift more severe in **lipid** cases
- SB plaque : Carina shift present only in the **lipid** case with **distal angle = 70°**
- **Lumen size of side branches marginally changes after stenting but lumen shape is effected**




Results: plaque shift

- Volumetric analysis




SB ostial compromise = lumen volume decrease in the SB ostial segment after MB stenting
Plaque shift = plaque volume increase in the SB ostial segment after MB stenting
Carina shift = vessel volume decrease in the SB ostial segment after MB stenting

- SB ostial compromise and plaque shift more severe in **lipid** cases
- Carina shift present only in the **lipid** case with **distal angle = 70°**
- **Lumen size of side branches marginally changes after stenting but lumen shape is effected**



Conclusions/Discussion (I)


- Development coronary bifurcation model including plaque dimensions
- Quantification of plaque shift by volumetric analysis → same order of magnitude as Xu et al. Circ cardiovasc int 2012
- Plaque shift depends on plaque composition and bifurcation angle
 - The largest plaque shift in lipid rich plaque
 - The largest plaque shift in small bifurcation angle




Conclusions/Discussion (I)

- Lumen size of SIDE BRANCH does marginally change --> more influence on lumen shape

Simulation



Xu et al. 2012



Acknowledgements



J.J. Wentzel, E. Regar, R. van Geuns, A. van der Steen, P. Serruys

Erasmus MC, Rotterdam, The Netherlands

C. Chiastra, F. Migliavacca, G. Dubini

Politecnico di Milano, Milano, Italy

F. Iannaccone, M. de Beule, P. Mortier

Gent University, Gent, Belgium




European Research Council

Established by the European Commission

The **10th** international symposium on
Biomechanics in Vascular Biology and Cardiovascular Disease



APRIL 23-24, 2015 ROTTERDAM THE NETHERLANDS



First Announcement



The
University
Of
Sheffield.

www.shearstresssymposium.nl